

1. Technical Data Sheet

Summary	TruCytes™ CD45(+), CD3(+), CD4(+), CD8(+), CD19(+), CD14(+), CD16(+), CD56(+) are lyophilized cell mimics to simulate TBNK phenotypes. They are formulated to work with lyse-wash, and lyse-no wash conditions.																																							
Application	This product is intended to provide positive signals for specified biomarkers in the table below:																																							
	<table border="1"> <thead> <tr> <th data-bbox="406 728 657 801">Antibodies</th> <th colspan="3" data-bbox="665 728 1415 801">Tested Clones</th> </tr> </thead> <tbody> <tr> <td data-bbox="406 801 657 864">CD45</td> <td data-bbox="665 801 916 864">2D1</td> <td data-bbox="924 801 1174 864">MEM-28</td> <td data-bbox="1182 801 1415 864">HI30</td> </tr> <tr> <td data-bbox="406 864 657 927">CD3</td> <td data-bbox="665 864 916 927">SK7</td> <td data-bbox="924 864 1174 927">*OKT3</td> <td data-bbox="1182 864 1415 927">UCH-T1</td> </tr> <tr> <td data-bbox="406 927 657 990">CD4</td> <td data-bbox="665 927 916 990">SK3</td> <td data-bbox="924 927 1174 990">OKT4</td> <td data-bbox="1182 927 1415 990">RPA-T4</td> </tr> <tr> <td data-bbox="406 990 657 1052">CD8</td> <td data-bbox="665 990 916 1052">SK1</td> <td data-bbox="924 990 1174 1052">HI8a</td> <td data-bbox="1182 990 1415 1052">RPA-T8</td> </tr> <tr> <td data-bbox="406 1052 657 1115">CD19</td> <td data-bbox="665 1052 916 1115">SJ35C1</td> <td data-bbox="924 1052 1174 1115">4G7</td> <td data-bbox="1182 1052 1415 1115">HIB19</td> </tr> <tr> <td data-bbox="406 1115 657 1178">CD16</td> <td data-bbox="665 1115 916 1178">B73.1</td> <td data-bbox="924 1115 1174 1178">3G8</td> <td data-bbox="1182 1115 1415 1178">CB16</td> </tr> <tr> <td data-bbox="406 1178 657 1240">CD56</td> <td data-bbox="665 1178 916 1240">NCAM-16</td> <td data-bbox="924 1178 1174 1240">MEM-188</td> <td data-bbox="1182 1178 1415 1240">MY31</td> </tr> <tr> <td data-bbox="406 1240 657 1296">CD14</td> <td data-bbox="665 1240 916 1296">M5E2</td> <td data-bbox="924 1240 1174 1296">61D3</td> <td data-bbox="1182 1240 1415 1296">HCD14</td> </tr> </tbody> </table>				Antibodies	Tested Clones			CD45	2D1	MEM-28	HI30	CD3	SK7	*OKT3	UCH-T1	CD4	SK3	OKT4	RPA-T4	CD8	SK1	HI8a	RPA-T8	CD19	SJ35C1	4G7	HIB19	CD16	B73.1	3G8	CB16	CD56	NCAM-16	MEM-188	MY31	CD14	M5E2	61D3	HCD14
Antibodies	Tested Clones																																							
CD45	2D1	MEM-28	HI30																																					
CD3	SK7	*OKT3	UCH-T1																																					
CD4	SK3	OKT4	RPA-T4																																					
CD8	SK1	HI8a	RPA-T8																																					
CD19	SJ35C1	4G7	HIB19																																					
CD16	B73.1	3G8	CB16																																					
CD56	NCAM-16	MEM-188	MY31																																					
CD14	M5E2	61D3	HCD14																																					
	* Very low to no signal. We do not recommend using this clone for this product																																							
	NOTE: Antibody clones not listed above need to be empirically tested to determine compatibility. It is highly recommended to titrate all antibodies in order to achieve the best results.																																							
	For Research Use Only. Not for use in diagnostic or therapeutic procedures.																																							
Materials	This product is lyophilized for stability and ease of use. Each vial contains approximately 2.5x10 ⁵ cell mimics.																																							
Handling and Safety	No special handling or safety precautions are necessary. See Safety Data Sheet (SDS) at www.slingshotbio.com .																																							
Storage	Store lyophilized products at -20°C upon receipt. Use immediately upon reconstitution.																																							
Expiration	One year from the date of manufacturing.																																							

Instructions for Use

1. Tap down the vial to ensure that all cell mimics are collected at the bottom of the vial.
2. Add 250 µL of 1x PBS buffer to the vial. Make sure to avoid contacting or disturbing the pellet until it has been rehydrated. Gently pipette up and down to resuspend the cell mimics. Transfer the content to an Eppendorf or FACS tube (other vessels such as plates can be used as well, and the volume can be adjusted accordingly).
3. Add 1mL of PBS to the original glass vial, mix, and transfer the remaining of the mimics to the tube prepared in the previous step for a final volume of 1.25 mL. Another 1 mL of 1x PBS can be added at this stage if the vessel is a FACS tube for a final volume of 2.25 mL. Centrifuge at appropriate speed (16000 x g for 2 minutes or 5000 x g for 5 minutes), remove the supernatant without disturbing the cell pellet.

Note: Washing with 2 mL of 1X PBS before staining is recommended to reduce non-specific binding.

4. Add an appropriate amount of your staining antibody cocktail and mix well by vortexing.

Note: Titrate antibodies on cells prior to making antibody cocktail for best results.

5. Incubate the mixture at RT in the dark for 10-15 min.
6. Add 2 mL of staining buffer for FACS tubes, or 200 µL if staining in wells (adjust according to the utilized vessel), vortex or mix, and spin down at 5000 x g for 5 minutes. If the staining is performed parallel to cells, the samples can be centrifuged at 400 x g for 5 minutes. Aspirate the supernatant without disturbing the pellet.

Note: For best results, stain in FACS tubes to ensure sufficient washing

7. Repeat the previous step x1.

Note: For best results, two to three washes are recommended.

8. Add desired volume of running or staining buffer to the FACS tube (or other utilized vessels).
9. Acquire cell mimics using the same instrument settings as leukocytes.

QC Data

Figure 1. Scatter plot and gating. The following shows the gating strategy for proper detection of subset populations. Image resolution is not looking good

